

ROTARY EVAPORATOR PILOT SCALE



Good Partner - Good Success

Ablaze Export manufacturer of large capacity rotary evaporator **GMP and Non GMP** designed for industrial plants Scale. It meets all the process needs by using Borosilicate glass and PTFE permitting aggressive substances treatment.

The unit works under vacuum and also under atmosphere conditions. The rotation of the **10 or 200 Liter Evaporation** flask in the thermostatic bath results in a forced convection and homogenous giving out of the product, thus preventing sedimentation.

Moreover, it ensures high evaporating exchange surface. Possibility to customize glassware as for clients needs .

❖ SALIENT FEATURES :

- ✓ Lift motorized system for heating Bath .
- ✓ Handling trolley for large capacity of flask .
- ✓ Stainless steel frame with wheels.
- ✓ GMP & NON GMP designed .
- ✓ High evaporation rate under vacuum.



Ablaze Export Incorporation
E-52/4, Sardar Estate , Ajwa Road
Vadodara 390019- Gujarat INDIA
www.ablazeexport.com